



Sterile, alcohol-based ready to use solution for disinfecting surfaces in cleanrooms.

perform® sterile alcohol EP

Our Plus

- fast, effective and residue free
- sterile (X-Ray irradiated) and microbial filtrated (0.2 µm)
- suitable for cleanrooms grade A/B, ISO class 5/6
- ready to use in a handy, sterile 500 ml spray system
- colour and perfume free
- efficacy tested according to European norms (EN)

Application areas

Disinfection of surfaces and equipment/devices in cleanrooms. Especially for Grade A and B clean areas, according to:

- EU Guide to Good Manufacturing Practice for medicinal products Annex 1: Manufacture of sterile medicinal products, March 2009.
- FDA Draft Guidance for Industry: Sterile Drug Products produced by aseptic processing - current Good Manufacturing practice, U.S. Department of Health and Human Services; Food and Drug Administration, September 2004.

Instructions for use

Apply perform® sterile alcohol EP undiluted to objects and surfaces and allow to take effect. Ensure complete wetting. The product dries off quickly. perform® sterile alcohol EP can be used for glove disinfection. Maximum application rate 50 ml/sq.m. Do not use to treat surfaces that are particularly sensitive to alcohol, such as acrylic glass.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal EN1276, EN13697	ready-to-use	1 min.

Efficacy	Concentration	Contact time
fungicidal EN1650, EN13697	ready-to-use	2 min.
virucidal against enveloped viruses in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	30 sec.



perform[®] sterile alcohol EP

Product data

Composition:
100 g solution contains the following active substances: 25 g Ethanol (94 %), 35 g Propan-1-ol.

Chemical-physical data

Color	colourless
Density	approx. 0,89 g/cm ³ / 20 °C
Flash point	27 °C / Method : DIN 51755 Part 1
Form	liquid
pH	approx. 7 / 20 °C

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.
For individual questions:
Application Department
Phone: +49 40 52100-666
E-Mail: info@schuelke.com

Special advice

Use disinfectants safely. Always read the label and product information before use.

The sterile 500 ml spray system offers the following advantages:

- convenient size for storage
- ergonomically shaped spray pump for ease of use
- remains sterile after opening

The bottles are sterile (with irradiation indicator) and doublewrapped individually. perform[®] sterile alcohol EP comprises alcohol of pharmaceutical quality in accordance with the Ph. Eur.

Information for order

Item	Delivery form	Item no.
perform [®] sterile alcohol EP 1 l bottle	6/Carton	on request
perform [®] sterile alcohol EP 500 ml bottle	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
perform [®] spray bottle holder	693540
perform [®] spray bottle holder - plate	693541

Related Products

- perform[®] sterile dry wipes
- perform[®] sterile mix dry wipes
- perform[®] sterile premium dry wipes

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com