

schülke -†



Wound irrigation solution for cleansing, moistening and decontamination of chronic skin wounds

octenilin® wound irrigation solution

Our plus:

- cleans wound areas that are difficult to access
- colourless and odourless
- excellent cleaning performance
- especially well tolerated by skin and tissue
- non irritating and non-sensitising
- for repeated and long term-use

Application areas

- rapid, effective cleansing and decontamination of wounds
- for the removal of wound crusts consisting of necrotic tissue, pathogens, biofilm and fibrinous films
- suitable for difficult to access areas, such as fissures and wound pockets
- gentle moistening of wounds
- moistening of dressings and wound pads
- gentle removal of difficult to remove, caked dressings/wound pads
- for combined use with Negative Pressure Wound Therapy (V.A.C. Ulta)

Application

Thoroughly irrigate and clean the wound with octenilin® Wound Irrigation Solution during each dressing change (octenilin® Wound Irrigation Solution can be warmed to body temperature immediately before application). Dressings and wound coverings soaked with octenilin® Wound Irrigation Solution can be used to dissolve fibrin coverings. The application should be repeated until all coverings and necroses can be removed and the wound is optically clean. Rinsing is generally not necessary. The wound can be further treated with octenilin® Wound Gel, if necessary.

Best before 8 weeks after first use.

Product profile

• Moisturising property

Due to its low surface tension, octenilin® Wound Irrigation Solution has excellent moisturising and cleansing properties on tissue surfaces, including difficult to reach areas, such as fissures and wound pockets.

• Cleansing power to assist removal of biofilms

In comparison to wound irrigation solutions such as Ringer or isotonic saline solutions commonly used in hospitals, octenilin® Wound Irrigation Solution shows significantly better cleansing power and, therefore, a better reduction of pathogens.

• Conservation

octenilin® Wound Irrigation Solution considerably reduces the risk of germs spreading to adjacent areas during application. If applied in accordance to instructions, octenilin® Wound Irrigation Solution offers an effective germ barrier between the wound and the patient.

• Tissue tolerability and biocompatibility

Non-irritating, non-sensitising, pain free, no tissue toxicity, no impairment granulation or epithelialisation of tissue. The good tolerability of octenidine has been confirmed by clinical experience over many years and by clinical studies on chronic wounds.

octenilin® wound irrigation solution

Product data

Ingredients:

Purified water, glycerol, ethylhexylglycerin, octenidin HCl.

Original condition:

Sterile, sealed. [STERILE/A](#)

Stability at room temperature, see stability information.

Use before:

8 weeks after first opening. The bottle must be closed after each use.

Chemical / physical data:

Appearance: colourless solution
Odour: odourless solution

Special advice

• Use with caked dressings and dressings that are difficult to remove

If the wound dressing is difficult to remove, soaking it with octenilin® Wound Irrigation Solution is recommended to facilitate the gentle removal of the dressing without traumatising the wound surface.

• Side effects

None observed.

• Contraindications

As allergic reactions can never be fully ruled out, octenilin® Wound Irrigation Solution should not be used if allergies to any of its contents are known or suspected.

octenilin® Wound Irrigation Solution should not be used on hyaline cartilage, the eyes, ears or nose, urinary bladder and in the abdominal cavity. To prevent possible tissue injury, the product must not be injected or applied to tissues pressure. Adequate drainage from wound cavities must be provided (eg. flexible drain tube).

• General safety instructions

- Only to be applied to the wound externally. Ensure that the preparation is not introduced under pressure into the tissue.
- Do not use for infusion or injection.
- Do not ingest.
- Only flawless and undamaged containers must be used.
- Do not expose to sunlight.
- Keep out of reach of children.

Information for order

Item	Delivery form	Item no.
octenilin® wound irrigation solution 350 ml	10/Carton	121 701

- ARTG 277 960

Related products

- octenilin® Wound Gel
- octenisan® Wash Lotion
- octenicare® Repair Cream

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.