



rotasept®

High performance, ready to use cleaning and disinfection solution developed for burs and other rotating instruments

Our Plus

- excellent material compatibility
- corrosion-preventing
- effective against bacteria (inc. TB) and viruses (inc. HBV, HCV, HIV and enveloped viruses) in 15 mins
- for quicker results insert the rotasept® bath into an ultrasonic bath
- ready to use

Application areas

For manual cleaning and disinfection of burs, cutters, diamond grinders and other rotating precision instruments. Due to its excellent material compatibility and comprehensive efficacy, rotasept® provides all important product properties for reprocessing rotating precision instruments as well as instruments made of steel and stainless steel.

Instructions for use

For professional use only. rotasept® is a ready-to-use product and is added undiluted to the bur bath.

If rotasept® is used as cleaner, place each small rotating instrument in the bur bath solution immediately after use in order to avoid surface-drying of dirt. Rinse all instruments thoroughly after reprocessing, using water of at least drinking water quality or preferably sterile distilled water or fully deionised water, to completely remove residues of the solution.

Immersion time for final disinfection: 15 minutes. Longer contact times should be avoided to ensure best material compatibility. Replace working solutions every working day and if contamination is clearly visible. National regulations may require that cleaning and final disinfection are carried out in two separate process steps. Do not mix with cleaning agents and disinfectants. Please follow the reprocessing recommendations of the instrument manufacturers.

Material compatibility:

Alkali- and alcohol-sensitive plastics, aluminium parts and instruments which have been made using synthetic resin adhesives should not be placed in rotasept®. We recommend gigasept® instru AF for these instruments.

Ultrasonic bath:

rotasept® is highly suitable for the use in ultrasonic bath.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>in accordance with VAH</i>	ready-to-use	15 min.
bactericidal EN13727, EN14561, <i>in accordance with VAH</i> - dirty conditions	ready-to-use	5 min.
yeastcidal EN13624, EN14562, <i>in accordance with VAH</i> - dirty conditions	ready-to-use	5 min.
yeastcidal - dirty conditions	ready-to-use	5 min.
virucidal against enveloped viruses <i>in accordance with DVV</i>	ready-to-use	5 min.
Vaccinia virus	ready-to-use	5 min.
virucidal <i>in accordance with DVV</i> (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	5 min.
virucidal EN14476, EN17111 - dirty conditions	ready-to-use	1 min.
Adenovirus	ready-to-use	5 min.
Poliovirus <i>in accordance with DVV</i>	ready-to-use	5 min.
Polyoma SV40	ready-to-use	5 min.

CE 0297



Product data

Composition:

100 g concentrate contains the following active substances:
1.6 g Potassium hydroxide, 10.0 g Propan-2-ol, 0.1 g 2-Ethylhexanol.

labelling according to Regulation (EC)648/2004: < 5 % non-ionic surfactants.

Special advice

Always read the label and product information before use.

Any serious incidents related to the product are to be reported to the manufacturer and competent authority. Keep container tightly closed. Protect from frost, heat and sunlight. Store in the original container. Do not use the product after the expiry date.

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit www.schuelke.com.

Information for order

Item	Delivery form	Item no.
rotasept -INT- 2 l FL	5/Carton	104741

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com