

# schülke →



## aspirmatic®

Cost effective, multi purpose cleaning, disinfection and de-odorising of dental suction systems and spittoons

### Our Plus

- non foaming solution
- fresh fragrance
- economical: makes 100L in-use solution
- simple and safe handling
- for amalgam traps of all kinds
- proven by extensive practical tests with the leading dental units

### Application areas

For daily cleaning, disinfection and de-odorising of the dental suction systems and spittoon.

### Product Profile

aspirmatic® (cleaning and disinfection) ensures extensive cleaning and the necessary hygienic safety during practice operation.

### Instructions for use

Use once or twice daily. For best results use aspirmatic® twice daily at the end of each clinical session. To make up the in-use solution add 40 ml aspirmatic® to 1960 ml luke warm water. To ensure the solution does not foam when diluting always measure the correct amount of luke warm water first, then pour the correct amount of aspirmatic® in to the water. Suck 1.5 litres of use solution through the suction hoses and pour the remaining 1/2 litre in the rinse basin. Leave to soak for minimum of 30 minutes, or overnight. After soaking always rinse through with at least one litre of clean water before using the suction system.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> <i>EN13727, EN14561</i>	2 % (20 ml/l)	30 min.
<b>bactericidal</b> <i>in accordance with VAH</i>	2 % (20 ml/l)	30 min.
<b>yeastcidal</b> <i>in accordance with VAH</i>	2 % (20 ml/l)	30 min.
<b>fungicidal</b> <i>EN13624, EN14562</i>	2 % (20 ml/l)	30 min.
<b>virucidal against enveloped viruses</b> <i>in accordance with DVV</i>	2 % (20 ml/l)	30 min.
<b>HBV</b>	2 % (20 ml/l)	30 min.
<b>HCV</b>	2 % (20 ml/l)	30 min.
<b>HIV</b>	2 % (20 ml/l)	30 min.
<b>Vaccinia virus</b>	2 % (20 ml/l)	30 min.

### Certificates

- VAH certificate

CE 0297



## Product data

### Composition:

100 g concentrate contains the following active substances:

7.5 g dioctyl dimethyl ammonium chloride.

Labelling according to Regulation (EC) No 648/2004: < 5 %

non-ionic surfactants, < 5 % soap, perfumes.

Other ingredients: corrosion inhibitors.

## Special advice

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit [www.schuelke.com](http://www.schuelke.com)

The UVV (accident prevention regulations) require the persons dealing with disinfectants which are not determined for the skin to wear gloves. Do not use after the expiration date.

## Information for order

Item	Delivery form	Item no.
aspirmatic -INT/NL- 2 I FL	5/Carton	150106

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail [mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



**schülke UK**  
Schülke & Mayr UK Ltd  
Cygnet House, Meadowhall  
Sheffield, S9 1AT  
Phone +44-114-254 35 00  
Fax +44-114-254 35 01  
[mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)