



For fast, effective and gentle wound irrigation

octenilin® wound irrigation solution

Our Plus

- excellent cleaning performance
- removes wound coating and biofilm
- moistens wounds and creates an ideal environment for wound healing
- preserves moisture of dressings and wound coverings
- especially well tolerated by the skin and tissues
- painless application
- suitable for long term use
- use for up to 8 weeks after first opening

Application areas

- rapid, effective cleansing of wounds
- for the removal of wound crusts consisting of necrotic tissue, pathogens, biofilm and fibrinous films
- suitable for difficult to access areas, such as fissures and wound pockets
- for moistening of wounds
- for moistening of dressings and wound pads
- for gentle removal of difficult to remove, caked dressings/ wound pads
- for combined use with Negative Pressure Wound Therapy (V.A.C.Ult, KCI)

Product Profile

Moistening property

Due to its low surface tension, octenilin® wound irrigation solution has excellent moistening properties and cleansing properties on tissue surfaces. The solution also works well in difficult to reach areas such as fissures and wound pockets.

Cleansing power to remove biofilms

Studies on the cleansing power against microorganisms in the biofilm revealed a good performance.

In comparison to other wound irrigation solutions, such as Ringer or isotonic saline solutions commonly used in hospitals, octenilin® wound irrigation solution shows a significantly better cleansing power and, therefore, a better reduction of pathogens.

Performance

octenilin® wound irrigation solution contains Octenidine as a preservative, capable of inactivating germs and thereby reducing the risk of germ transmission into the surrounding area during application. octenilin® wound irrigation solution offers an effective germ barrier between the wound and the patient if applied correctly.

Tissue tolerability and biocompatibility

Non-irritating, non-sensitising, painless, without tissue toxicity, no impairment granulation and no impairment of the epithelialisation. The good tolerability of Octenidine has been confirmed over many years clinical experience and clinical studies on chronic wounds.

Instructions for use

Irrigate and clean the wound with octenilin® wound irrigation solution during each change of dressing (octenilin® wound irrigation solution can be warmed up to body temperature immediately before application). Dressings can also be removed this way. Dressings and wound coverings soaked with octenilin® wound irrigation solution can be used to dissolve fibrin coatings. The application should be repeated until all coatings and necroses can be removed and the wound is visibly clean. Rinsing afterwards is generally not necessary. The wound can be further treated with octenilin® wound gel, if necessary.

CE 0297



octenilin® wound irrigation solution

Product data

Ingredients:

Aqua valde purificata, Glycerol, Ethylhexylglycerin, Octenidine HCl

Original condition:

Sterile, sealed. Stability at room temperature, see stability information.

Use before:

8 weeks after first opening. The bottle must be closed after each use.

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit www.schuelke.com

Use with caked dressings, i. e. dressings that are difficult to remove.

If the wound dressing is difficult to remove, soaking it with octenilin® Wound Irrigation Solution is recommended to facilitate the gentle removal of the dressing without traumatising the wound surface.

Side effects

None observed.

Contraindications

octenilin® Wound Irrigation Solution should not be used in patients with known or suspected allergy to one or more of the ingredients. If in doubt consult a physician. octenilin® wound irrigation solution should not be used on hyaline cartilage, the eyes, ears, nose, urinary bladder and in the abdominal cavity! To prevent possible tissue injury, the product must not be injected or applied to tissues with pressure. Adequate drainage from wound cavities must be provided (e.g. for flexible drain tube).

General safety instructions

Do not use for infusion or injection. Not for oral use. Only to be applied to the wound externally. Use only intact and undamaged packs. Protect the product against exposure to direct sunlight. To prevent the introduction of bacteria when using octenilin® wound irrigation solution, ensure that the packaging does not come into contact with the wound. Medical device - Keep out of the reach and sight of children.

Order reference codes:

Cat. code: MRB443

PIP code: 348-9069

Information for order

Item	Delivery form	Item no.
octenilin wound irrigation 350 ml FL	10/Carton	121701

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnets House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com