



Powdered concentrate based on active oxygen for disinfecting surfaces.

perform® classic concentrate OXY

Our Plus

- sporizidal and material compatible
- convenient single-use sachets
- colour free
- efficacy tested according to European norms (EN)

Application areas

Disinfection of all types of surfaces. perform® classic concentrate OXY is suitable for use in all high-risk areas where, in addition to broad effectiveness, strict requirements are placed on the environmental behaviour and user-friendliness of a surface disinfectant. Due to the excellent material compatibility of the active solution, perform® classic concentrate OXY can also be used with problematic materials, such as acrylic glass (incubators).

Instructions for use

To prepare the working solution, first add water (preferably hand-warm) to the tub or bucket, then sprinkle in and dissolve the product.

For example:

1 sachet (40 g) for 4 L = 1 % solution

Heavily contaminated surfaces are to be pre-cleaned. Ensure complete wetting. Contact time has to be considered. Rinse surfaces intended to be in direct contact with final products. To prevent the build-up of residue on treated surfaces, keep residual moisture on surfaces to a minimum after wiping. perform® classic concentrate OXY has excellent cleaning properties, therefore mixing it with cleaning boosters is not necessary. The composition of the active ingredient system in the active solution remains stable for one working day. If the active solution becomes very dirty, replace it earlier.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>EN1276, EN13697</i>	0.5 % (5 g/l)	30 min.
fungicidal <i>EN1650, EN13697</i>	1.5 % (15 g/l)	60 min.
virucidal against enveloped viruses <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	0.5 % (5 g/l)	5 min.
Vaccinia virus	0.5 % (5 g/l)	5 min.
sporizidal <i>EN13697</i>	1 % (10 g/l)	2 h
sporizidal <i>EN13704</i>	2 % (20 g/l)	2 h



perform[®] classic concentrate OXY

Product data

Composition:

100 g powder contains the following active substances: 45 g Pentapotassium bis(peroximonosulphate) bis(sulphate).
Labelling according to Regulation (EC) No. 648/2004: 5 - 15 % anionic surfactants, < 5 % non-ionic surfactants, < 5 % soap, < 5 % phosphonates, perfumes.

Chemical-physical data

Color	white
Flash point	Not applicable
Form	granular
pH	ca. 4 / 5 g/l / 20 °C / in water
Viscosity, dynamic	Not applicable

Special advice

Use disinfectants safely. Always read the label and product information before use.

Do not store bags above room temperature. Wearing of protective gloves (e.g. made of Nitrile rubber) is recommended.

The active solution may bleach textile dyes that are highly sensitive.

Information for order

Item	Delivery form	Item no.
perform [®] classic concentrate OXY 250 x 40 g	250/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com