



Sterile ready to use solution based on peracetic acid for disinfecting surfaces in cleanrooms.

## perform® sterile PAA

### Our Plus

- fast, effective and sporicidal
- sterile (γ-Ray irradiated)
- ready for use in the handy 500 ml spray system or in the 300 ml bottle with screw cap for use in combination with the perform® sterile dry wipes
- colour and perfume free
- efficacy tested according to European norms (EN)

### Application areas

Disinfection of surfaces and equipment / devices in cleanrooms, in particular for class A and B cleanroom areas, according to:

- Annex 1 of EU Guide to Good Manufacturing Practice: Manufacture of sterile medicinal products, March 2009.
- FDA Guidance for Industry: Sterile drug products produced by aseptic processing - current Good Manufacturing Practice; U.S. Department of Health and Human Services, Food and Drug Administration, September 2004.

### Instructions for use

- Apply perform® sterile PAA 500 ml spray bottle undiluted to pre-cleaned objects and surfaces and allow to take effect, ensuring that all surfaces are moistened completely.
- perform® sterile PAA 300 ml bottle with screw cap for use in combination with perform® sterile dry wipes. Ensure complete wetting and keep the surface moist during the entire application time. Ensure that all visible soiling is removed before disinfection.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> EN1276, EN13697	ready-to-use	5 min.

Efficacy	Concentration	Contact time
<b>fungicidal</b> EN1650, EN13697	ready-to-use	5 min.
<b>virucidal against enveloped viruses</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	1 min.
<b>sporicidal</b> EN13697, EN13704	ready-to-use	5 min.



## Product data

Composition:  
100 g solution contains the following active substances: 0.09 g peracetic acid.

### Chemical-physical data

Color	colourless
Density	1,01 g/cm <sup>3</sup> / 20 °C
Flash point	> 105 °C
Form	liquid
pH	ca. 3,5 / 10 g/l / 20 °C / in water
Viscosity, dynamic	1,05 mPa*s / 20 °C

## Special advice

**Use disinfectants safely. Always read the label and product information before use.**

Professional use only. Read safety data sheet before use.

Good material compatibility with stainless steel, PTFE, PVDF, PEHD, PP, Viton, glass, etc.

Contact with non-ferrous metals (brass, copper) should be avoided. Wearing of protective gloves (e.g. made of Nitrile rubber) is recommended. Carry and store product container in an upright position.

## Information for order

Item	Delivery form	Item no.
perform <sup>®</sup> sterile PAA 500 ml bottle	10/Carton	on request
perform <sup>®</sup> sterile PAA 300 ml bottle	10/Carton	on request
perform <sup>®</sup> sterile PAA 300 ml bottle	10/Carton	on request
perform ster PAA nonirr 300 ml FL	10/Carton	on request
perform sterile PAA -n.i. EU- 500 ml FL	10/Carton	on request
perform sterile PAA -n.i. CN- 500 ml FL	10/Carton	on request
perform sterile PAA -EU- 500 ml FL	10/Carton	on request
perform sterile PAA -EU- 500 ml FL	10/Carton	on request
perform sterile PAA -CN- 500 ml FL	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Application aids

Application aids	Item no.
perform <sup>®</sup> spray bottle holder	693540
perform <sup>®</sup> spray bottle holder - plate	693541

## Related Products

- perform<sup>®</sup> sterile dry wipes
- perform<sup>®</sup> sterile mix dry wipes
- perform<sup>®</sup> sterile premium dry wipes

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)

Manufacturer:

Bioxal

ZI Sud Secteur A

Route des Varennes

71100 Chalon sur Saône

France



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

Manufacturer  
Bioxal  
ZI Sud Secteur A  
Route des Varennes  
71100 Chalon sur Saône  
France

schülke Headquarters  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)