



dosit® des

Dispenser for the production and application of aqueous disinfectant solutions.

Optionally a second or third alternative concentration can be adjusted

Our plus:

- Microprocessor-controlled system
- Monitoring of complete dispensing procedure by sensors
- Dispensing range: 0.25 % – 5 % in 0.05 % increments adjustable; dispensing tolerance + 5 %
- Solution discharge depressurised directly at pivoting outlet
- Automatic switch-off at lack of concentrate or water as well as at disruption of concentrate flow
- Clear text display of all operating conditions incl. the adjusted concentrations
- Mains supply isolation according to DIN 1988, Part 4/TRWI) / EN 1717
- Built according to RKI guidelines RKI test under preparation

Description

dosit® des is a microprocessor-controlled unit for the production of ready-to-use disinfectant up to 5 %. All conventional disinfectants with a viscosity of < 50 mP*s can be used.

- The unit is designed for decentralized, stationary applications and discharges the solution for use depressurised.
- The operating conditions are shown on a 32-digit clear text display.
- Operation and programming takes place via a membrane keypad.
- The disinfectant concentrates can be adjusted in 0.05 % increments.
- The solution submission is between 1l - 5l in 0.5 l steps; between 5l - 50l in 1 l steps. The solution submission > 50l occurs with shutting down automatically


- The housing is made of stainless steel in a stylish bevelled finish. The suction lance and canister bracket can be mounted on the right or left of the unit. The canister console is easy to remove for cleaning purposes.
- The DGWV tested mains supply isolation enables its direct connection to a drinking water connection.
- The unit is built and tested according to all current regulations (VDE [German Association for Electrical, Electronic and Information Technologies], EMC [Electromagnetic compatibility], DVGW [German Association for Gas and Water], Robert Koch Institute [RKI], RKI test under preparation).

Technical Data


Technical Data	
Dispensing quantity:	7 l/min
Dispensing range:	0.25 % – 5 %
Dispensing tolerance:	+ 5 %
Cold water connection:	3/8" male thread, brass
Water inlet pressure:	1.5 – 10 bar
Power supply:	230 V/~ 50 Hz
Output:	max. 75 Wat
Connecting cable:	1.5 m with safety plug
Dimensions:	
H	578 mm
B	510 mm
T	180 mm
Tare weight:	12 kg

Description

- | | |
|---|--------------------|
| 1 | Dispensing pump |
| 2 | Suction lance |
| 3 | Pressure regulator |
| 4 | Dirt trap |
| 5 | Water inlet |
| 6 | Solution outlet |
| 7 | Dispensing nipple |
| 8 | Control system |



Produkt	Article no.
dosit [®] des	644000



Mounting instructions

Electrical connection 1.5 m standard
equipped with top Schuko plug.
Left or right mounted.

Optional accessories	Article no.
Canister holder lockable dosit [®] des	644017
Console for 10-l-canister	6000055



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com