



Exclusively for the testing of glove permeability.

## schülke optics lotion

### Our Plus

- fluorescent solution for varification of glove permeability
- shows the permability of gloves by means of UV-light

### Application areas

Exclusively for the testing of glove permeability.  
Do not use as a care lotion.

### Product Profile

Fluorescent solution only for varification of glove permeability.

### Instructions for use

Put the gloves on. Rub schülke optics lotion into gloved hands in the same way as for hygienic or surgical hand disinfection. Remove the gloves and hold the hands under a UV Lamp (366 nm). Areas of skin which were affected by the lotion will glow. This indicates gaps in the gloves.



# schülke optics lotion

## Product data

Aqua, Paraffinum Liquidum, Cetearyl Ethylhexanoate, Glycerin, Simmondsia Chinensis (Jojoba) Seed Oil, Polyglyceryl-3 Methylglucose Distearate, Stearic Acid, Glyceryl Stearate, Isopropyl Myristate, Cera Alba, Allantoin, Sodium Polyacrylate, Xanthan Gum, Stilbene Derivative, Potassium Sorbate, Sodium Benzoate, Citric Acid

## Chemical-physical data

Color	white
Density	approx. 0,98 g/cm <sup>3</sup> / 20 °C
Flash point	Not applicable
Form	emulsion
pH	approx. 5,0 / 20 °C
Viscosity, dynamic	approx. 1.000 mPa*s / 20 °C / Method : DIN 53019

## Information for order

Item	Delivery form	Item no.
schülke optics lotion 500 ml FL	20/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)