



Low-alcohol, rapid disinfection of non-invasive medical devices and surfaces.

## mikroqid® universal liquid

### Our Plus

- suitable for disinfecting sensitive communication devices
- dermatologically tested
- ready to use

### Application areas

The product is suitable for disinfecting medical devices and other surfaces. Our versatile, rapid disinfectant is also suitable for disinfecting sensitive communication devices. Smartphones and tablets are among the most common hand contact surfaces. With mikroqid® universal liquid, these can be safely prepared without damaging the materials.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> EN13727 - dirty conditions	ready-to-use	15 sec.
<b>bactericidal</b> EN16615 - clean conditions	ready-to-use	1 min.
<b>bactericidal</b> in accordance with VAH - clean conditions	ready-to-use	1 min.
<b>tuberculocidal</b> EN14348 - clean conditions	ready-to-use	5 min.
<b>yeastcidal</b> in accordance with VAH - dirty conditions	ready-to-use	2 min.
<b>yeastcidal</b> EN16615 - clean conditions	ready-to-use	1 min.
<b>yeastcidal</b> in accordance with VAH - clean conditions	ready-to-use	2 min.
<b>yeastcidal</b> EN13624 - dirty conditions	ready-to-use	1 min.

Efficacy	Concentration	Contact time
<b>virucidal against enveloped viruses</b> EN14476 - dirty conditions	ready-to-use	15 sec.
<b>limited spectrum virucidal activity</b> EN14476 - dirty conditions	ready-to-use	15 min.
<b>Norovirus</b> EN14476 - dirty conditions	ready-to-use	1 min.
<b>Norovirus</b> EN14476 - clean conditions	ready-to-use	30 sec.
<b>Polyoma SV40</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	15 sec.
<b>Rotavirus</b> EN14476 - clean conditions	ready-to-use	15 sec.

### Certificates

- VAH certificate
- IHO listed



## Product data

100 g solution contains the following active substances: 17,4 g Propan-2-ol, 12,47 g Ethanol (94 % w/w)  
 Labelling according to Regulation (EC) No. 648/2004: < 5% anionic surfactants

### Chemical-physical data

Color	colourless
Density	ca. 0,95 g/cm <sup>3</sup>
Flash point	26 °C / Method : DIN 51755 Part 1
Form	liquid
pH	3 - 3,6
Viscosity, dynamic	not determined

## Special advice

**Use disinfectants safely. Always read the label and product information before use.**

The product has very good material compatibility with metals and plastics. Keep away from heat and sources of ignition. Further information is available on request. The shelf-life of the opened container can be found on the packaging.

### Instruction for use

Apply undiluted to the surface, wipe and allow taking effect. Ensure complete wetting and maintain the surface moist during the entire exposure time. Ensure that all visible dirt is removed before disinfection. Spray the surface from a distance of approx. 30 cm (approx. 40 ml/m<sup>2</sup>) and allow drying. For adenovirus and Mycobacterium tuberculosis, maintain the surface wet during the entire exposure time. **For professional use only.** Any serious incidents related to the product are to be reported to the manufacturer and competent authority. Not for final disinfection of semicritical and critical medical devices! The product can be used without gloves, provided that infection prevention and occupational safety allow this.

## Information for order

Item	Delivery form	Item no.
mikrozid® universal liquid 1 l spray bottle (RFL)	10/Carton	on request
mikrozid® universal liquid 1 l bottle (WSP)	10/Carton	on request
mikrozid universal liquid 750 ml FL	1/Carton	on request
mikrozid® universal liquid 5 l	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810
Decanter funnel	117901
FOAM Spray pump DIN 32	180123
measuring cup 500 ml	136101
Opener all-in-one (1l eurobottle/can/container)	70000804
schülke dosing feeder 5 l / 10 l (20ml / stroke)	117101
tap for 5 l / 10 l can	135502

## Related Products

- mikrozid® AF liquid
- mikrozid® sensitive liquid
- mikrozid® universal wipes premium

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:  
 Application Department  
 Phone: +49 40 52100-666  
 E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
 Schülke & Mayr GmbH  
 Robert-Koch-Str. 2  
 22851 Norderstedt, Germany  
 Phone +49 40 - 52100 - 0  
 Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)