



Powdered concentrate for disinfecting surfaces / PT2

perform® classic concentrate PEROXY

Our Plus

- sporicidal and very good material compatibility
- convenient single-use sachets
- colour and perfume free
- suitable for cleanroom classes C/D and non-classified areas
- efficacy tested according to European norms (EN)

Application areas

For the disinfection of surfaces, perform® classic concentrate PEROXY is suitable for use in all risk areas, where, in addition to broad efficacy, high demands are placed on the environmental behavior and the user-friendliness of a surface disinfection preparation. Due to the very good material compatibility of the working solution perform® classic concentrate PEROXY can also be used on problematic materials, such as acrylic glass (incubators), can be used.

Instructions for use

To prepare the working solution, first add water (max. 30 °C) to the tub or bucket, then sprinkle in perform® classic concentrate PEROXY and dissolve. For example: 1 sachet (40 g) to 4 liters of water = 1 % (10 g/l). Remove visible soiling beforehand. Wipe surfaces well to wet them and observe the reaction time. If necessary, rinse the surface in contact with the product. To reduce the formation of residues on treated surfaces, surfaces should only be left with minimal residual moisture after wiping. perform® classic concentrate PEROXY has good cleaning properties. Mixing with cleaning boosters is therefore not allowed. The prepared working solutions have a lifetime of max. 30 hours (at 20 °C) and 4 hours (at 40 °C) for guaranteed stability of the active ingredient. The working solution should be replaced earlier in case of heavy contamination.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal EN1276, EN13697	0.5 % (5 ml/l)	30 min.
yeastcidal EN1650, EN13697	0.5 % (5 ml/l)	60 min.
fungicidal EN1650, EN13697	1.5 % (15 ml/l)	60 min.
virucidal against enveloped viruses EN14476, EN16777	0.5 % (5 ml/l)	5 min.
virucidal EN14476, EN16777	1 % (10 ml/l)	5 min.
Bacillus subtilis EN13704	1 % (10 ml/l)	120 min.



perform[®] classic concentrate PEROXY

Product data

100g contains: 49.5 g pentakalium bis (peroxymonosulfate) bis (sulfate).

Chemical-physical data

Colour: white

Flash point: not applicable

Form: solid powder

pH value: ca. 3,0 / 20 g/l / 20 °C / in Wasser

Viscosity, dynamic: not applicable

Special advice

Use disinfectants safely. Always read the label and product information before use. For professional use only.

The wearing of protective gloves (e.g. made of nitrile rubber) is recommended. Do not store bags above room temperature.

The active solution may bleach textile dyes that are highly sensitive.

Information for order

Item	Delivery form	Item no.
perform [®] classic concentrate PEROXY 40 g	60/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- perform[®] sterile concentrate PEROXY

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com