

schülke -†



3L decontamination bath

bath-system 3l

Our Plus

- branded bath with inner instrument tray and lid



bath-system 3l

Product data

Internal dimensions: 250 x 150 x 120 mm
External dimensions: 300 x 185 x 125 mm

Special advice

Max. temperature 55 °C/ 131 °F

Information for order

Item	Delivery form	Item no.
set - bath 3 lt - cover white SE	1/Carton	144307

Environmental information

schulke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com>. For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.