



Small decontamination bath for burs and other rotating instruments

## bur bath

### Our Plus

- use with rotasept® to clean and disinfect burs and other rotating instruments
- can be used as a static bath or inserted into an ultrasonic bath
- 250 ml capacity
- non-autoclavable



## Information for order

Item	Delivery form	Item no.
schülke bur box	12/Carton	144030

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com>. For any queries you might have about the product or using it please call 01 14 2543500 or E-mail [mail.uk@schuelke.com](mailto:mail.uk@schuelke.com).



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.