



Pink nitrile examination glove made from an exceptionally soft material for maximum comfort

Sempercare[®] Nitrile Soft Gloves

Our Plus

- textured finger tips for excellent grip even when moist
- extremely comfortable to wear
- safe alternative for type I allergy sufferers
- pink coloured glove

Product Profile

Not only is the innovative, powder-free nitrile examination glove incredibly soft, the pink colour also brings colour to the treatment room. As well as providing protection from infections, this comfortable glove offers excellent tactility and is suitable for users and patients with type I allergies. Compliant with all European impermeability standards for proven protection against liquids and viruses. Textured finger tips guarantee a sure grip even in a moist environment.

CE 0121



Sempercare[®] Nitrile Soft Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves made of a nitrile butadiene rubber (NBR), powder-free, non-sterile. Ambidextrous. Free from natural latex Synthetic copolymer. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations. For single use. Check for damage before use, do not use damaged gloves.

Information for order

| Item | Delivery form | Item no. |
|---|---------------|----------|
| Soft Pink PFN Gloves - XL 180 ST BX | 10/Carton | 352074 |
| Soft Pink PFN Gloves - Medium 200 ST BX | 10/Carton | 352072 |
| Soft Pink PFN Gloves - Large 200 ST BX | 10/Carton | 352073 |
| Soft Pink PFN Gloves - Small 200 ST BX | 10/Carton | 352071 |
| Soft Pink PFN Gloves - XSmall 200 ST BX | 10/Carton | 352070 |

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com