



High quality soft nitrile gloves with easy donning material

Protects Clinic Nitrile Gloves *Super Light*

Our Plus

- powder-free, non sterile gloves
- safe for allergy type I sufferers

Product Profile

Examination gloves made of a nitrile, powder -free, non-sterile material. The gloves contain no latex proteins making them a safe alternative for type I allergy sufferers. The ambidextrous design ensures that the gloves are suitable for all hand types, and they are compliant with Dir. 93/42/EEC, EN 455.



Protects Clinic Nitrile Gloves

Super Light

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

This product contains dithiocarbamate and zinc mercaptobenzothiazol - not to be used in the case of hypersensitivity to these substances. Formulation available on request. For single use. Check for damage before use, do not use damaged gloves. Do not use for chemical hazards. Store in original packaging in a dry and dark places at 10° to 30°C. Dispose of in accordance with valid regulations.

Information for order

Item	Delivery form	Item no.
protects superlight PFN Large 200 ST BX	10/Carton	70000152
protects superlight PFN XL 180 ST BX	10/Carton	70000153
protects superlight PFN Small 200 ST BX	10/Carton	70000150
protects superlight PFN XS 200 ST BX	10/Carton	70000149
protects superlight PFN Med 200 ST BX	10/Carton	70000151

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



EMAS

GEPRÜFTES
UMWELTMANAGEMENT
REG. NR. 01-100-10003



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnets House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com