



Powder free, non-sterile examination glove with extended cuff

Sempercure[®] Nitrile Shine+ Gloves

Our Plus

- extra long cuff (265mm) offers extra safety to wrists and forearms
- increased material wall thickness increases glove strength and user safety
- innovative foil opening for quick and easy glove removal

Product Profile

The nitrile examination glove offers an extra-long cuff for additional user safety by covering wrists and lower arms. This powder and latex free glove has a higher wall thickness, and is a safe alternative for latex type I sufferers. 150 appealing neutral white colour gloves in each box.

CE 0121



Sempercare[®] Nitrile Shine+ Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves made of a nitrile butadiene rubber (NBR), powder-free, non-sterile. Ambidextrous. Free from natural latex Synthetic copolymer. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations. For single use. Check for damage before use, do not use damaged gloves.

Information for order

Item	Delivery form	Item no.
Nitrile Shine+ Gloves Medium 150 ST BX	10/Carton	70000408
Nitrile Shine+ Gloves Small 150 ST BX	10/Carton	70000407
Nitrile Shine+ Gloves Large 150 ST BX	10/Carton	70000409
Nitrile Shine+ Gloves XL 150 ST BX	10/Carton	70000410
Nitrile Shine+ Gloves Xsmall 150 ST BX	10/Carton	70000406

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnets House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com