



Comfortable high quality latex gloves with easy donning

Sempercare® Latex Gloves

Our Plus

- soft latex material to fit the shape of your hand perfectly
- micro-rough surface for excellent touch and safe grip
- inner coated for improved donning
- natural white latex colour

Product Profile

Textured finger tips for good handling of instruments and an excellent sense of touch. The latex material is very soft and elastic giving the wearer a high level of comfort. The gloves have an ambidextrous shape and rolled cuffs to help ease donning.

CE 0121



Sempercare[®] Latex Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves from made latex proteins, powder-free, non-sterile. Ambidextrous. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations.

This product contains latex proteins - not to be used in the case of hypersensitivity to these substances. For single use. Check for damage before use, do not use damaged gloves.

Information for order

Item	Delivery form	Item no.
Sempercare PFL Gloves - Med 100 ST BX	10/Carton	352035
Sempercare PFL Gloves - XL 90 ST BX	10/Carton	352036
Sempercare PFL Gloves - Small 100 ST BX	10/Carton	352033
Sempercare PFL Gloves - Large 100 ST BX	10/Carton	352034

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



EMAS
GEPRÜFTES
UMWELTMANAGEMENT
REG. NR. D-150-10003



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com

 **Air Liquide** A company of the Air Liquide-
HEALTHCARE Group