



High quality soft nitrile gloves with easy donning

Sempercure® Nitrile Gloves

Our Plus

- good resistance to chemical substances and punctures
- ambidextrous with rolled cuffs for ease of fitting
- suitable for type I allergy sufferers
- proven protection against harmful microorganisms
- blue coloured glove

Product Profile

These powder-free gloves do not contain any natural latex proteins. This makes them perfect for users with type I allergies who do not want to miss the greater comfort of elastic natural latex. The blue colour makes it easier to distinguish from latex gloves. The soft material increases wearing comfort and makes the gloves comfortable and safe to work with, even for longer periods of time. The structured finger tips guarantee a good grip, even in moist conditions. The higher puncture resistance and improved chemical resistance with the same thickness as natural latex, offers the user protection and the feeling of safety.

CE 0121



Sempercare[®] Nitrile Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves made of a nitrile butadiene rubber (NBR), powder-free, non sterile. Ambidextrous. Compliant with Dir. 93/42/EEC, EN 455, PPE 89/686/EEC, EN 420, EN 374 (deviating length because examination glove as well 240+/-5mm). Free from natural latex Synthetic copolymer. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations.

This product contains dithiocarbamate and zinc mercaptobenzothiazolate - not to be used in the case of hypersensitivity to these substances. For single use. Check for damage before use, do not use damaged gloves.

Information for order

Item	Delivery form	Item no.
Sempercare PFN Gloves - XL 90 ST BX	10/Carton	352040
Sempercare PFN Gloves - Large 100 ST BX	10/Carton	352039
Sempercare PFN Gloves - Small 100 ST BX	10/Carton	352037
Sempercare PFN Gloves - Med 100 ST BX	10/Carton	352038

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



EMAS

GEPRÜFTES
UMWELTMANAGEMENT
REG. NR. 0-100-00023



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com

 **Air Liquide** A company of the Air Liquide-
HEALTHCARE Group