



Soft nitrile examination glove combine excellent sensitivity and high elasticity

Sempercare® Nitrile Velvet Gloves

Our Plus

- very high elasticity
- extra soft material, feels like a second skin
- excellent sense of touch and superior wearing comfort
- safe alternative for latex type I allergy sufferers
- lavender blue colour

Product Profile

Sempercare® Velvet is powder-free and does not contain any natural rubber protein. As such, this glove is perfectly suitable for glove users with Type-I allergy.

The innovative material formulation comes with a reduced wall thickness which allows for maximum protection against infections while providing superior comfort. The velvety lavender-blue material wears like a second skin and is especially gentle on sensitive skin.

CE 0121



Sempercare[®] Nitrile Velvet Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves made of a nitrile butadiene rubber (NBR), powder-free, non-sterile. Ambidextrous. Free from natural rubber protein. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations. For single use. Check for damage before use, do not use damaged gloves.

Information for order

Item	Delivery form	Item no.
Sempercare velvet PFN - Small 200 ST BX	10/Carton	70001708
Sempercare velvet PFN - XS 200 ST BX	10/Carton	70001707
Sempercare velvet PFN - XL 180 ST BX	10/Carton	70001712
Sempercare velvet PFN - Large 200 ST BX	10/Carton	70001711
Sempercare velvet PFN - Med 200 ST BX	10/Carton	70001709

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke UK
Schülke & Mayr UK Ltd
Cygnets House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com