



Alcohol-based ready to use solution for disinfecting surfaces.

## perform® classic alcohol EP

### Our Plus

- fast, effective and residue free
- ready to use
- colour and perfume free
- efficacy tested according to European norms (EN)

### Application areas

For disinfection of surfaces in areas where fast microbiological efficacy is required.

### Instructions for use

Apply perform® classic alcohol EP undiluted to objects and surfaces and allow to take effect. Ensure complete wetting. The product dries off quickly. Maximum application rate 50 ml/sq.m. Do not treat surfaces that are particularly sensitive to alcohol such as acrylic glass.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> <i>EN1276, EN13697</i>	ready-to-use	1 min.
<b>fungicidal</b> <i>EN1650, EN13697</i>	ready-to-use	2 min.
<b>virucidal against enveloped viruses</b> <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	ready-to-use	30 sec.



# perform<sup>®</sup> classic alcohol EP

## Product data

### Composition:

100 g solution contains the following active substances: 25 g Ethanol (94 %), 35 g Propan-1-ol.

### Chemical-physical data

Color	colourless
Density	0,890 g/cm <sup>3</sup> / 20 °C
Flash point	27 °C / Method : DIN 51755 Part 1
Form	liquid
pH	approx. 7 / 20 °C

## Special advice

**Use disinfectants safely. Always read the label and product information before use.**

## Information for order

Item	Delivery form	Item no.
perform <sup>®</sup> classic alcohol EP 1 l bottle	10/Carton	on request
perform <sup>®</sup> classic alcohol EP 200 l FA	1/drum(s)	on request
perform <sup>®</sup> classic alcohol EP 10 l canister	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810
Compressed-air nozzle	120001
Decanter funnel	117901
schülke tap for 5 l / 10 l can	135501

## Related Products

- perform<sup>®</sup> classic mix dry wipes

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)