



For the application of hand disinfectants, washing and care lotions for 500 ml and 1.000 ml disposable bottles.

Preparation dispenser Type KHL

Our Plus

- meets the requirements of the trade association (BGR 250 / TRBA 250) and the workplace regulations for direct removal of products
- adjustable dosing stroke of approx. 0.75 ml to approx. 1.5 ml
- casing aluminium, anodised (bacteriostatic effect)
- autoclavable and dishwasher-safe dispenser
- easy assembly
- available as a counter-dispenser (incl. counting function)

Application areas

Easy and exact dosing of hand disinfecting, washing and care preparations.

- Dosing of the required preparations via a dosing block at the dosing pump to amounts of approx. 0.75 ml to 1.5 ml per stroke.
- The exchangeable pump is made of stainless steel and ensures a smooth preparation flow. No clogging, clotting or after-dripping.
- Casing made of matt silver anodised aluminium. Arm lever made of stainless steel.
- If required, the dispenser can be removed from the wall mounting plate, e.g. for cleaning.
- Quick filling with the original bottle.
- Easy withdrawal of the dosing pump via the front.



Preparation dispenser Type KHL

Information for order

Item	Delivery form	Item no.
Product dispenser long arm 500 ml	1/Carton	on request
dispenser KHL 500 - stainl. steel pump	1/Carton	on request
Präparatespender KHL 500 Counter	1/Carton	on request
Product dispenser long arm 1000 ml	1/Carton	on request
dispenser KHL 1000 - stainl. steel pump	1/Carton	on request
Präparatespender KHL 1000 Counter	1/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Schalenhalter f.1 L Spender	669390
Schalenhalter f.500ml Spender	669240

Related Products

- desderman® pure
- hyclick® System
- Preparation dispenser KHN touchless
- preparation dispenser sm2
- Preparation dispenser Type KHK



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com