

# schülke -t



Phosphoric acid-based neutralizer for automated reprocessing of medical equipment.

## mucapur<sup>®</sup> N

### Our Plus

- removes residual alkalinity from the cleaning step
- prevents tarnishing on metallic equipment, corrosion and scale deposits in machines and items for rinsing
- low residue behaviour (does not influence sensitive microbiological laboratory testing)
- free of tensides
- mucapur<sup>®</sup> N is recommended by Miele

### Application areas

- neutralisation of medical instruments, particularly metallic equipment, after alkaline cleaning
- suitable for the use in all washer/disinfectors and multi-chamber washer/disinfectors with liquid dosage systems

### Instructions for use

Concentration:  
Starting at 0.1 %. Dosing is performed by machine-integrated dosing pumps. Please note recommendation of machine manufacturer.

CE



## Product data

Components:  
anorganic acidic

### Chemical-physical data

Color	colourless
Density	1,33 g/cm <sup>3</sup> / 20 °C
Flash point	Not applicable
Form	liquid
pH	< 2 / 20 °C / concentrate
Viscosity, dynamic	Not applicable

## Special advice

Contains phosphoric acid.

## Information for order

Item	Delivery form	Item no.
mucapur N 200l	1/drum(s)	on request
mucapur N 30l	1/Canister	on request
mucapur-N 5 l KA	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Accessories

Accessories	Item no.
Can key for 5 + 10 l	135810

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com). For individual questions: Customer Care Phone: +49 40 52100-666  
E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



A company of the Air Liquide-Group

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 (0) 40 - 52100 - 0  
Fax +49 (0) 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[mail@schuelke.com](mailto:mail@schuelke.com)

Schülke & Mayr UK Ltd.  
1 Jenkin Road  
GB-Sheffield S9 1 AT  
Phone +44 (0) 1142 - 5435 - 00  
Fax +44 (0) 1142 - 5435 - 01  
[www.schuelke.co.uk](http://www.schuelke.co.uk)  
[mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

Schülke & Mayr Ges.m.b.H.  
Seidengasse 9  
1070 Wien, Austria  
Phone +43 (0) 1 - 5232501-0  
Fax +43 (0) 1 - 5232501-60  
[www.schuelke.at](http://www.schuelke.at)  
[office.austria@schuelke.com](mailto:office.austria@schuelke.com)