

# schülke -†



Tenside-based rinse aid and drying accelerator for automated reprocessing of medical equipment with neutralizing effect

## mucapur<sup>®</sup> S

### Our Plus

- improves speed of water drainage and shortens drying time by up to 60 % - which saves energy, time and money
- guarantees dryness after cleansing
- neutralizing effect removes any remaining alkalinity following cleansing
- provides incredible shine with no staining; cares and protects against deposits in the machine
- very economical due to low levels of concentration needed for effectiveness
- mucapur<sup>®</sup> S is recommended by Miele

### Application areas

- special clear rinsing agent for medical use in washer/disinfectors and multi-chamber washer/disinfectors

### Instructions for use

Concentration:  
starting at 0.01 % to 0,03 %. Dosing is performed by machine-integrated dosing pumps. Please note recommendation of machine manufacturer.

CE



## Product data

### Composition:

5 - 15 % nonionic surfactants, < 5 % anionic surfactants

Further ingredients: organic acids, adjuvants

### Chemical-physical data

Color	colourless
Density	1,05 g/cm <sup>3</sup> / 20 °C
Flash point	47,5 °C
Form	liquid
pH	< 3 / concentrate
Viscosity, dynamic	No data available

## Information for order

Item	Delivery form	Item no.
mucapur S 5000 ml	3/Carton	on request
mucapur-S 1l FL	10/Carton	on request
mucapur-S 200 l FA	1/drum(s)	on request
mucapur-S 5 l KA	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Accessories

Accessories	Item no.
Can key for 5 + 10 l	135810

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com). For individual questions: Customer Care Phone: +49 40 52100-666 E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



A company of the Air Liquide-Group

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 (0) 40 - 52100 - 0  
Fax +49 (0) 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[mail@schuelke.com](mailto:mail@schuelke.com)

Schülke & Mayr UK Ltd.  
1 Jenkin Road  
GB-Sheffield S9 1 AT  
Phone +44 (0) 1142 - 5435 - 00  
Fax +44 (0) 1142 - 5435 - 01  
[www.schuelke.co.uk](http://www.schuelke.co.uk)  
[mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

Schülke & Mayr Ges.m.b.H.  
Seidengasse 9  
1070 Wien, Austria  
Phone +43 (0) 1 - 5232501-0  
Fax +43 (0) 1 - 5232501-60  
[www.schuelke.at](http://www.schuelke.at)  
[office.austria@schuelke.com](mailto:office.austria@schuelke.com)