



For the determination of the total plate count.

mikrocount® TPC

Our Plus

- Easy to use
- No preparation time
- Testing of surfaces, semi-solid materials and liquids with just one product
- Expensive laboratory equipment is not necessary
- Results within 24 - 72 hours

Application areas

This nutrient media allows the growth of the most common bacteria, yeasts and fungi. The agar contains lecithin and Tween for neutralising any potential residues of disinfectants with cationic substances.



Chemical-physical data

Color	not determined
Density	Not applicable
Flash point	Not applicable
Form	Solid form
pH	Not applicable
Viscosity, dynamic	Not applicable

Information for order

Item	Delivery form	Item no.
mikrocount TPC 20 ST BX	18/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- mikrocount[®] duo
- mikrocount[®] TPC / E



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.