



Dosing container for preparing and pouring out aspirmatic® and aspirmatic® cleaner

## schülke dental application aid

### Our Plus

- a very good 2 liter scale helps to prepare the stock solution
- left over stock solution for the mouth rinse basin stays in the dosing container
- easy attachment of the extraction hose on to available adapters

### Application areas

The dosing container is for preparing and pouring out the stock solution into the extraction system. It makes working with aspirmatic® and aspirmatic® cleaner very easy and safe. As the scale is very good to see the preparation of the stock solution is made easy. As the rest of the stock solution stays automatically in the dosing container it can be used for the mouth rinse basin. The easy attachment of the extraction hoses due to the available adapters makes the individual working steps easy.



# schülke dental application aid

## Information for order

Item	Delivery form	Item no.
schülke dental application aid	1/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Related Products

- aspirmatic®
- aspirmatic® cleaner

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product:

[www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)